

GERSTEL

wister[®]
THE ORIGINAL



GERSTEL Twister®




The GERSTEL Twister® enables efficient extraction of organic compounds from aqueous matrices based on Stir Bar Sorptive Extraction (SBSE). SBSE is a solvent-free extraction technique, which is significantly faster than most conventional extraction techniques. SBSE is up to 1000 x more sensitive than SPME since the stir bar has significantly more sorbent volume and since it can extract, and concentrate analytes from, a much larger sample volume due to the efficient stirring. The Twister may look like a conventional magnetic stirring rod, but while it stirs samples, such as water, body fluids or beverages it absorbs and concentrates organic compounds into its sorbent coating. A large number of samples can be extracted simultaneously using multi-position stir plates resulting in high productivity and throughput. Labor- and resource-intensive sample preparation steps are eliminated. Analytes are typically desorbed from the Twister using thermal desorption. Liquid desorption using a solvent can be performed whenever the analytes are very high-boiling, thermally labile, or when they must be determined by HPLC.

The GERSTEL Twister® offers a simple possibility to position one or more Twisters magnetically on the inner wall of a sample vial for more efficient sample extraction. The Twister approach enables the extraction of a sample using more Twisters simultaneously, for example using Twisters with different types of phases for a more complete combined extraction and improved analyte recovery.



Twister/SBSE Benefits:

- Up to 1000 times more sensitive than SPME
- Quantitative, with large linear range
- A large number of samples can be extracted simultaneously for highest productivity
- Requires minimal time and labor
- Thermal desorption and GC/MS analysis are performed in one automated system



Heated stir plates used to extract 100 mL water samples with Sediment. 100 compounds are determined in one GC/MS run at ultra-low concentrations (EU-WFD).

Twister Sorbent Phases

	PDMS Twister *	EG/Silicone Twister
		
Phase	Polydimethylsiloxane (PDMS)	Polydimethylsiloxane (PDMS) / Ethylene glycol (EG) - copolymer on an inert metal grid for mechanical stabilization
Enrichment	Unspecific sorption of apolar compounds $\log(K_{ow}) > 4$ The polarity range can be extended by adding salt to the sample (salting out).	Unspecific sorption of apolar compounds, additionally specific binding of polar hydrogen bond donors, such as phenols
Application examples	Pesticides in water PAH in marine tissue 2,4,6-TCA in wine Flavor compounds in food	Flavor compounds in beverages

* Available in 10 mm or 20 mm length with phase thickness 0.5 or 1.0 mm. Phase volume: 24, 47, 63, and 126 μ L respectively.

Twister Application Fields

The Twister has proven itself in a wide range of application areas, such as food and beverages, flavor and fragrance, extractables and leachables from packaging as well as environmental analysis. Thanks to its mechanical stability, the Twister can be reused up to 200 times, depending on the sample matrix and desorption conditions, and can also be used for field sampling; it is easy and convenient to use, requires little or no solvent and it is inexpensive to ship to the laboratory for analysis.

Applications:

- Food and beverages
- Flavor and Fragrance
- Environmental analysis of water (EU-WFD) or waste water
- Extractables and leachables from pharmaceutical packaging



www.gerstel.com/en/apps-twister-sbse.htm

GERSTEL Flex-Twister



The Flex Twister was developed specifically for Solvent Assisted Stir Bar Sorptive Extraction (SA-SBSE), a novel extraction technique that relies on Flex Twisters that have been swollen with a user defined solvent.

The solvent is added to the Flex Twisters before the extraction step leading to a modification of the sorbent phase properties and thereby, depending on the solvent used, to significantly improved extraction recovery for polar compounds with $\log K_{ow}$ values between 1.0 and 2.0. Compared with standard SBSE, the SA-SBSE phase volume is increased significantly leading to improved recovery even of non-polar compounds as an added bonus.

Following SA-SBSE, the analyst can desorb the Flex Twister using either a solvent for subsequent GC/MS or LC/MS analysis or thermal desorption with moderate temperature programming allowing the solvent to be evaporated before analytes are transferred to the GC/MS for determination.

GERSTEL Twicester®



The GERSTEL Twicester® offers a simple possibility to position one or more Twisters magnetically on the inner wall of a sample vial for more efficient sample extraction. The Twicester approach enables the extraction of a sample using more Twisters simultaneously, for example using Twisters with different types of phases for a more complete combined extraction and improved analyte recovery. The method, referred to as Multi-Stir Bar Sorptive Extraction (^mSBSE), was developed in a cooperation between GERSTEL K.K., Tokyo, Japan, and the Research Institute for Chromatography (RIC) of Professor Pat Sandra in Kortrijk, Belgium. Twicester is well suited for the simple and efficient determination of a wide range of flavor compounds in aqueous matrices and it has been used successfully for the analysis of beverages. A key element of the method is the use of both PDMS- and EG-Silicone Twisters in one vial and the simultaneous desorption of both Twisters combined into a single GC/MS run. While the PDMS Twister stirs the sample extracting non-polar to medium polarity compounds, the EG-Silicone Twister is held in position on the inside wall of the vial by the Twicester clip.

The EG-Silicone Twister very efficiently extracts more polar compounds that form hydrogen bonds as hydrogen donors, for example phenols. Twicester can be used to place Twisters in the liquid phase or in the headspace as required.

GERSTEL Twicester® benefits

Magnetic mounting of the Twister on the inside wall of the vial


- The sorbent surface is stationary avoiding mechanical stress.
- No sample to sample carry-over, only the Twister is in contact with the sample.

Multi-Stir Bar Sorptive Extraction (^mSBSE)

- Simple solution for simultaneous extraction of a sample with multiple Twisters based on the GERSTEL Twicester®
- A wide polarity range can be covered by using different Twisters simultaneously.
- Improved recovery and lower limits of determination through use of multiple Twisters

Simultaneous thermal desorption of multiple Twisters

- All extracted compounds are determined in one comprehensive GC/MS run
- Improved limits of determination through simultaneous desorption and analysis of multiple Twisters



SA-SBSE is a rugged analysis technique that is well suited for analysis of samples with a heavy matrix load, such as, flavor compounds in smoothies (GERSTEL AppNote 198).



Magnetic positioning of up to three Twisters using Twicester



^mSBSE using two or more Twisters



Transfer of the Twisters to a TDU Tube



Simultaneous thermal desorption of the Twisters followed by Cryofocusing in the CIS, and GC/MS analysis

Tube Conditioner TC



Clean and well-conditioned TD tubes and Twisters are essential in order to fully take advantage of the high sensitivity of the thermal desorption technique and to achieve the best possible quality of analysis. The TC 2 operates independently of the GC/MS system, eliminating the risk of system contamination during conditioning and keeping the GC/MS available for analytical tasks ensuring maximum productivity. During conditioning in the TC 2, TD tubes are either held at a user-specified fixed temperature or undergo a specified temperature program cycle while being purged with inert gas to efficiently remove contaminants. Up to 5 Twisters will fit into an empty TDS tube, allowing up to 50 Twisters to be conditioned in one batch. Each tube has independent gas flow regulation enabling the conditioning of one to ten tubes without the need to adjust flows. The TC 2 helps to ensure that conditioned TD tubes and Twisters are always available providing best possible quality of analysis.



Twister Desorption in the TDU



The Thermal Desorption Unit (TDU 2) is a flexible automated solution for thermal desorption and thermal extraction. The TDU 2 fits on top of any modern GC without the need for additional bench space and it is perfectly suited for the analysis of gaseous, liquid and solid samples. The TDU 2 incorporates the latest advances in thermal desorption technology. Intelligently designed and based on a "Liner-in-Liner" concept it has no valves or transfer lines. The TDU 2 is connected directly to the GERSTEL Cooled Injection System (CIS), which serves both as a cryo-focusing trap and as a temperature programmable GC inlet. Active sites are eliminated, reducing the risk of analyte loss, discrimination and memory effects to an absolute minimum. The TDU 2 can be operated in split or true splitless mode enabling it to cover the widest range of concentrations, to protect the column from water and contamination and to achieve the lowest possible detection limits. The TDU 2 low-flow split pneumatics provide improved flexibility and performance. For extreme sensitivity, the multi-desorption mode can be selected by mouse-click in MAESTRO.

Twister Back Extraction TBE



Twister® Back Extraction as MPS option

The combination of Twister technology (SBSE; Stir Bar Sorptive Extraction) with solvent elution and subsequent chromatographic determination is a powerful and reproducible process for the determination of high boiling and thermally labile organic analytes from aqueous samples; e.g. PAH, pesticides, herbicides or phenols.

The process was developed by UFZ Leipzig-Halle GmbH in close cooperation with GERSTEL.

Twister Back Extraction TBE offers:

- Extraction, enrichment and analysis of high boiling or thermally labile compounds
- Combination of Stir Bar Sorptive Extraction (SBSE) and LC or LC/MS analysis

GERSTEL MAESTRO Software

Next generation software for automated sample preparation and sample introduction. MAESTRO optimizes performance and throughput of GERSTEL modules and systems.

- Stand-Alone operation, fully integrated in the Agilent ChemStation and MassHunter Software, or integrated with the Thermo Scientific® Xcalibur™ sequence table
- Sample Prep by Mouse-Click using PrepBuilder functions
- Scheduler for easy planning of sequences and of laboratory work-flow
- PrepAhead / Multi-Sample Overlap: Automated overlapping of sample preparation and analysis for maximum throughput
- Priority samples can be added to the system at any point in the analysis sequence
- LOG file and Service LOG file functions ensure traceability
- Automated E-mail notification if the sequence is stopped
- Real-time monitoring of all modules and parameters
- Interactive on-line help function



www.gerstel.com

GERSTEL

G L O B A L A N A L Y T I C A L S O L U T I O N S

GERSTEL, Inc., USA
+1 410 - 247 5885
sales@gerstelus.com

GERSTEL BRASIL
+55 11 5665 8931
gerstel_brasil@gerstel.com

GERSTEL GmbH & Co. KG,
Germany
+49 208 - 7 65 03-0
gerstel@gerstel.com

GERSTEL AG, Switzerland
+41 41 - 9 21 97 23
gerstelag@ch.gerstel.com

GERSTEL K.K., Japan
+81 3 57 31 53 21
info@gerstel.co.jp

GERSTEL Shanghai Co. Ltd
+86 21 50 93 30 57
china@gerstel.com

GERSTEL LLP, Singapore
+65 6779 0933
sea@gerstel.com



Agilent Technologies
Premier Solution Partner