
**GERSTEL**

SamplePrep-  
Solution  
**3-MCPD**

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[1] Lucas, D.; Hoffmann, A.; Gil, C. Fully Automated Determination of 3-MCPD and Glycidol in Edible Oils by GC/MS Based on the Commonly Used Methods ISO 18363-1, AOCS Cd 29c-13, and DGF C-VI 18 (10) GERSTEL AppNote No. 191, 2017

GERSTEL Application Note No. 191, 2017

**3-MCPD and Glycidol contamination.** Table 1 shows the results from assay B, listing the amount of 3-MCPD determined in three different edible oil samples as well as the reference values.

**Table 1. 3-MCPD amount found in three different edible oils in mg/kg.**

3-MCPD	Amount [mg/kg]	
	Reference	Automated
Oil 1	0.77	0.68
Oil 2	0.68	0.63
Oil 3	0.27	0.29

For a given edible oil sample, the difference between the results for assays A and B multiplied by the previously determined conversion factor is used to calculate the amount of Glycidol in the sample. In table 2, the amounts obtained using this method are listed along with reference values.

**Table 2. Glycidol amount found in three different edible oils in mg/kg.**

Glycidol	Amount [mg/kg]	
	Reference	Automated
Oil 1	0.14	0.12
Oil 2	0.44	0.31
Oil 3	0.11	0.06

To demonstrate the good repeatability of the automated sample preparation method, five samples of the same edible oil were analyzed undergoing individual sample preparation and analysis. Table 3 shows the repeatability based on the entire sample preparation procedure and the subsequent GC/MS analysis.

**Table 3. Repeatability for 3-MCPD and Glycidol (n=5 samples).**

#	Amount [mg/kg]	
	3-MCPD	Glycidol
1	0.77	0.14
2	0.68	0.12
3	0.68	0.44
4	0.27	0.11
5	0.29	0.06

**CONCLUSIONS**  
In this work, we have shown that method ISO 18363-1 can be automated using the GERSTEL MPS and that the results obtained correlate well with reference data. This method is similar to two other frequently used methods: AOCS Cd 29c-13 and DGS C-VI 18 (10). The excellent relative standard deviations achieved for the complete process including GC/MS analysis speak in favor of the presented automation solution.

The work presented here involves an automated evaporation step as prescribed in the abovementioned official methods. This ensures that for most matrices, the required limits of detection can be reached using a single quadrupole mass spectrometer (MSD). A further important aspect of the evaporation step is that it removes excess derivatization reagent, which could otherwise build up in the GC/MS system and influence system stability.

**OUTLOOK**  
The described automation steps are not limited to the presented method. Such methods have already been tested for derivatization methods like the recently presented 3 in 1 approach, and can be adapted for that method with similar performance. The presented method has the advantage of being able to analyze a sample for Glycidol, 3-Monochloropropanediol (3-MCPD) and additionally 2-Monochloropropanediol (2-MCPD), all in a single run.

In addition to extracting and determining MCPD and Glycidol esters, the described automation platform can also extract and determine PAHs from edible oils using automated solid phase extraction combined with GC/MS determination.

**LITERATURE**  
[1] <http://www.bfr.bund.de>

GERSTEL Application Note No. 191, 2017

**RESULTS AND DISCUSSION**

**Chromatogram m/z 198:** Top: Virgin olive oil used as blank oil. Middle: Edible oil sample B (3-MCPD). Bottom: Edible oil sample assay A (3-MCPD + Glycidol).

The linearity of the method was verified by analyzing virgin olive oil spiked at five different levels. This was performed for both assays. In figure 5, the excellent linearity ( $R^2 > 0.9998$ ) achieved for both assays from 0.12 – 1.9 mg/kg is shown.

**Conditions:**  
3 µL injection volume  
baffled liner, deactivated solvent vent  
40°C (0 min); 12°C/s; 300°C (5 min)  
30 m Rxi-17 sil ms (Restek)  
d = 0.25 mm d<sub>i</sub> = 0.25 µm  
He, constant flow = 1 mL/min  
50°C (2 min); 10°C/min; 200°C (20°C/min; 300°C (5 min)  
Selected ion monitoring SIM  
3-MCPD: 196/198/147 and 3-MCPD-d5: 201/203/150 and

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GERSTEL®, GRAPHPACK®, Twicester and TWISTER® sind eingetragene Marken der GERSTEL GmbH & Co. KG. Änderungen vorbehalten. Copyright by GERSTEL GmbH & Co. KG
**GERSTEL, Inc., USA**
**GERSTEL GmbH & Co. KG, Deutschland**
**GERSTEL AG, Schweiz**
**GERSTEL Co. Ltd, Shanghai**